

MA-2 Micromanipulator Stand

For holding the manipulator.


This stand is for temporary storage of stereotaxic manipulators when not in use. Specially designed to accommodate the manipulators' particular shape.



* <u>SM-15R</u> sold separately.

AP frame bar size

Narishige produces world standard 18.7mm square shape AP frame bar. Please be advised when mounting on your existent frames since our stereotaxic manipulators are tailored for 18.7mm AP frame bar.





MA-2 Micromanipulator Stand

Specification

Dimensions/Weight	W250 × D150 × H85mm, 2kg

Specifications are subject to change without notice. Copyright © NARISHIGE Group. All Rights Reserved. 1999-2019